


Knipsa

Endomyocardial Biopsy Forceps



- Sharp edges
- No squeeze artifacts
- Allows the collection of tissue samples

Single use Biopsy Forceps

Further advantages

Excellent maneuverability

- easy to navigate
- Bowden cable is easy to move

Quick and efficient

- ambulant and stationary
- low complication rate



Order information

| REF | Usable length | Cutting jaws size | Cutting jaws volume | Recommended sheath size | Packaging |
|-----------|---------------|-------------------|----------------------|-------------------------|------------------------|
| Knipsa500 | 50 cm | 1.80 mm | 2.42 mm ³ | 6 Fr (2.0 mm) | sterile, single packed |

Knipsa is for single use only. The Biopsy Forceps from Joline is available in a standard size and used in cardiology.

Joline GmbH & Co. KG
Neue Rottenburger Str. 50
72379 Hechingen
Germany

phone: +49 (0) 7471 9881-0
fax: +49 (0) 7471 9881-222
e-mail: info@joline.de
www.joline.de



Quality first: Joline is DIN EN ISO 13485:2016 certified.
All Joline products are CE-Marked.
Subject to alterations: 51499/07/20_RV003